

# GeneExpressMap Quickstart Guide

## *Getting Started*

### **Requirements**

GeneExpressMap is a MATLAB GUI application and as such requires a working installation of MATLAB

### **Installation**

Click the "Download Source" button above and extract the archive into a folder of your choosing. Start MATLAB and change its working directory to the folder in which you have extracted the downloaded archive. At the MATLAB prompt, type "GeneExpressMap" to begin!

## *Example Run*

### **Loading data sets**

On opening, GeneExpressMap displays four empty image windows on the right side of the GUI as well a set of processing buttons and settings on the left side. GeneExpressMap is designed to operate on three sets of data per experiment. The first data set is a nuclear stain, the second and third are fluorescent whole mount *in situ* hybridization (F-WMISH) data sets. To populate the GUI with data, select **Start New Analysis** from the **File** menu and load tif stacks of the appropriate data when prompted. If you wish to analyze a single channel F-WMISH data set, simply select the F-WMISH to be analyzed data twice.

### **Setting Analysis Parameters**

Once the nucleus and FISH data sets are loaded into GeneExpressMap, the default settings for all analysis parameters are set for the GUI. The user must set the **X/Y Resolution** and **Z Resolution** for their data set. These values will be unique to each dataset depending on how the data was acquired. For **Threshold Level**, these settings are calculated for each data set. For all other parameters, the defaults are general settings that have yielded satisfactory results on a variety of data sets. The user may start with these settings and refine them later if desired.

### **Image Processing**

A typical run through GeneExpressMap will next proceed with image thresholding, denoising, and dilation of the data sets in order to produce an overlap between the nuclear stain and F-WMISH signals. These steps are carried about by pressing the threshold, denoise, and dilate buttons as they become active on the left hand side of the GUI.

### **Nucleus identification and F-WMISH Expression assignment**

A model of all nuclei in the data set is next constructed by pressing the **Find Nuclei** button. This constructs a three dimensional model of all nuclei from the the two dimensional stack of processed images in the nuclear stain data set.